

Installing WinISI / WinScan updates.

*****Warning. These files are for non-translated versions of WinScan or WinISI. Installation over translated versions of these program could result in overwriting of the translated text files. If you have translated copies, please contact your local FOSS product specialist or click Feedback on the website for further information.***

Please follow the instructions below to update your WinISI or WinScan installation:

1. Back up your existing WinISI or WinScan installation. The ISI software is usually located in C:\winisi or C:\winscan. Copy this entire directory to a safe place on your hard drive or a network drive.
2. Download the update file to the desktop or other directory on the computer hard drive.
3. Use WinZip to open the Zip file, and then double click the .exe file to automatically extract the contents of the ZIP file to your existing WinISI or WinScan directory (typically c:\winisi). In response to the prompt, answer Yes to overwrite the old files.

The WinISI / WinScan installation is now updated.